

Introducing a Breakthrough in Skin Protection for Amputees...



The GlideWear® Prosthetic Liner Patch provides a protective low friction interface between amputees' residual limb and prosthetic liner; delivering immediate relief and lasting protection from discomfort, redness, and skin breakdown.

Indications

The GlideWear Prosthetic Liner Patch is indicated for amputees who experience skin irritation and areas of shear-related skin breakdown.

This product is not indicated for placement directly against open wounds. Consult your primary caregiver for wound dressing recommendations to accompany use of the GlideWear Prosthetic Liner Patch.

Features & Benefits

Easy to Use

Place GlideWear® Prosthetic Liner Patch directly against the area of concern, before donning liner over the top of the patch and limb.

Use a clean patch daily to protect against future breakdown.

Stays in Place all Day Long

GlideWear® Prosthetic Liner Patches are held in place by the gel liner or prosthetic socket; ensuring uninterrupted skin protection all day long.

Reduces Shear

Patented dual-layer GlideWear® fabric technology absorbs friction and related shear forces at the skin - prosthetic liner interface.

Durable & Reusable

Machine or hand wash and dry frequently to keep clean and dry.

Specifications

Construction: 85% Nylon, 15% Spandex.

Available in Two Sizes: Small (2.5 in. x 4.75 in.) *Part no. GW-PRO-2SB*
Large (4.5 in. x 7.5 in.) *Part no. GW-PRO-2LB*

For more information, visit www.GlideWear.com/PLP



Distributed Worldwide by:

